

HEXAGON HIV

Certified Quality – Recognised Worldwide

- > Excellent sensitivity
- > Meets high European standards
- > Fast, simple testing procedure
- > Storage at room temperature
- > Long shelf life

RAPID SCREENING TESTS

CoreLab DX



Rapid Screening Tests

Human

Diagnostics Worldwide

HEXAGON HIV

Rapid Screening Test for Whole Blood, Serum or Plasma Application



HEXAGON HIV

- > 3rd generation rapid test for the detection of antibodies to HIV
- > HIV-1 and HIV-2 antibodies are indicated with separate test lines

EXCELLENT PERFORMANCE

- > 100% sensitivity with HIV positive specimens (700 specimens including all subtypes)
- > 99.9% specificity with 1,552 HIV negative samples covering healthy individuals and hospital patients
- > State of the art sensitivity with HIV-1 serconversion panels



MEETS HIGH EUROPEAN STANDARDS

- > Additional assurance of quality – a benefit for users worldwide
- > Constant high quality – all new lots are tested by an independent external laboratory prior to release

MEETS MAJOR REQUIREMENTS

- > Suited for use with serum, plasma and whole blood
- > Storage at 2...30 °C – no need to maintain a cold chain
- > 24-month shelf life from the date of manufacture

		Cat. No.
HEXAGON HIV	40 tests	57002P
	100 tests	57004P

FAST AND SIMPLE TESTING PROCEDURE

