

HumaNex A1c

Automated HPLC HbA1c Determination

- > HbA1c analysis in less than 3 minutes
- > 50 samples onboard capacity
- > NGSP certified, CV < 3%

Diabetes Care



Human


Diagnostics Worldwide

HumaNex A1c

Gold Standard Technology for HbA1c Testing

Fully automated system

- > HbA1c analysis in less than 3 minutes
- > Throughput of 22 tests per hour
- > Onboard capacity of 50 samples per autoloader
- > Automatic sample identification
- > Automatic mixing
- > Sample volume:
Diluent mode 10 µl
Whole blood mode 1.5 ml



Intuitive user interface

- > LCD display
- > Simple touchscreen operation


Performance characteristics

- > Direct determination of HbA1c
- > Imprecision: CV < 3 %
- > Measuring range: 4 % to 17 %
- > Reporting in IFCC and DCCT values
- > NGSP certified

Ordering information

	REF
HumaNex A1c	16180
Reagent kit (800 tests)	16185
Calibrator kit	16187
Control kit	16189

Excellent correlation



Accuracy: HumaNex A1c correlates excellent with reference HPLC system

HumaNex A1c system reagents

